



Able Spacer with

Positive Whistle Mask

Technical specifications

Tradition and Innovation - Visionary thinking and a fascination with technology have guided us throughout our history to develop innovative products of outstanding reliability. Anticipating trends to improve the quality of life.

 **CLEMENT CLARKE
INTERNATIONAL**

Able Spacer with Positive Whistle Mask

Technical specifications



Positive Whistle Mask		(with Able Spacer)
Dimensions	73 x 86 x 42mm	73 x 86 x 197mm
Weight	21g	83.1g

Product	Part Codes
Able Spacer with small mask (boxed)	3609121*
Whistle Mask	(tba)

* UK pack includes Whistle Mask

Product information

Able Spacer (Small Mask variant) is now available with the New Positive Whistle Mask, this guides drug delivery and works with the Rafi-Tone Lite App. The whistle in the mask will not sound if it does not fit the face properly, which makes it easy to confirm the fit. The Rafi-Tone Lite App can be downloaded free from Apple App Store and Google Play Store. The Rafi-Tone App works with the whistle tone as an incentive tool for children to have fun with, while learning to accept the use of their spacer. Each time the whistle sounds the game moves on.

The Able Spacer is a popular valved inhalation holding chamber that facilitates the use of pMDI medication. This improved version of the Able Spacer is now manufactured in the UK, made from transparent plastic incorporating anti-microbial silver ion technology which provides improved protection from contamination, and inhibits microbial growth. The visible valve movement makes it easier for parents or carers to observe and confirm correct inhalation and actuation.

Positive Whistle Mask

- Confirms correct fit
- Encourages correct flow rate

Positive Whistle Mask / Rafi-Tone App

- Encourages child acceptance of face mask
- Rafi-Tone Lite App is available as a free download

Able Spacer

- Accepts most pMDIs available
- Inhibits microbial growth
- Low-static properties